


# Dental Handpiece Lubrication Chart

Lubricate after each use and before sterilization. Never use disinfectant or disinfectant wipes on handpieces as they are corrosive!

## HIGHSPEED HANDPIECES

**Do Not Oil Front End**




### **For Aerosol Application:**

1. Fit spray can with appropriate adapter nozzle for handpiece.
2. With blank bur installed, spray into drive air line for 1 second (see illustration below). Hold handpiece and nozzle together while spraying, pointing front of head away from your face.
3. Attach handpiece to air supply. Cup head with paper towel, being careful not to touch bur. Run handpiece for 20 seconds, expelling lubricant/cleaner onto towel. Repeat if expelled lubricant is discolored.
4. Remove bur, clean handpiece with alcohol wipe and sterilize.

### **For Dropper Application:**

1. With blank bur installed, apply 3 drops into drive air line (see illustration below).
2. Attach handpiece to air supply. Cup head with paper towel, being careful not to touch bur. Run handpiece for 20 seconds, expelling lubricant onto towel. Repeat if expelled lubricant is discolored.
3. Remove bur, clean handpiece with alcohol wipe and sterilize.



## ATTACHMENTS & REMOVABLE HEADS

1. Use dropper to lubricate all moving parts & run briefly.
2. Remove bur, clean handpiece with alcohol wipe and sterilize.

## SLOWSPEED HANDPIECES

### **For Aerosol Application:**

1. Fit spray can with appropriate adapter nozzle for handpiece.
2. With blank bur installed, spray into drive air line for 1 second (see illustration). Hold handpiece and nozzle together while spraying, pointing front of head away from your face.
3. Attach handpiece to air supply. Cup head with paper towel, being careful not to touch bur. Run handpiece for 20 seconds, expelling lubricant/cleaner onto towel. Repeat if expelled lubricant is discolored.
4. Remove bur, clean handpiece with alcohol wipe and sterilize.

### **For Dropper Application:**

1. Apply 2 drops into the drive air line (see illustration). Run for 30 seconds.
2. Remove bur, clean handpiece with alcohol wipe and sterilize

## ELECTRIC HANDPIECES

1. Fit Spray can with appropriate adapter nozzle for handpiece.
2. Remove handpiece motor and bur from handpiece.
3. Slide lubrication adapter into back of handpiece. Spray for 2 seconds, holding handpiece and nozzle together and pointing front of head away from your face.
4. Attach handpiece motor, replace the bur and run for 20 seconds to purge. Repeat if expelled lubricant is discolored.
5. Remove bur, clean handpiece with alcohol wipe, and sterilize.